

2/28/2006

Elisha Wood-Charlson

Rodriguez-Lanetty, M., Krupp, D.A. and Weis, V.M. (2004) *Mar Ecol Prog Ser* 275: 97-102.

Hawaii – Infecting *Fungia scutaria* with Algae (quantitative measure of infection success)

Materials/Reagents

0.25µm filtered seawater (fsw)

Hemocytometer, coverslips, cell counter

6-well culture dish

Frozen *Artemia* sp.

40µm cell strainers (BD Falcon: 352340, blue)

1.5ml tubes

1ml plastic pipettes

10X PBS

Make up in 800 ml	800 ml
0.02M NaH ₂ PO ₄ (monobasic) 120 g/m	1.92
0.077M Na ₂ HPO ₄ (dibasic) 142 g/m	8.8
1.4 M NaCl 58 g/m	69.6
add DH ₂ O to 800 ml	
Autoclave	

1X PBS

100 ml 10X PBS

900 ml DH₂O water

Paraformaldehyde fix

1 g paraformaldehyde

12.5 ml DH₂O

2.5 ml 10X PBS

10 ml DH₂O

- Heat 12.5 DH₂O in a small beaker with vigorous stirring.
- Weigh out paraformaldehyde making sure to wear a mask and use gloves.
- Add paraformaldehyde to fast stirring, heated water – stir for about a minute, do not allow to boil. Powder will not dissolve.
- Remove from heat and continue to stir.
- Add 1 drop of 1 N NaOH to solution. The solution should clear.
- Once clear, filter the solution through a sterile filter and check that volume is 12.5 ml – if not add up to 12.5 with water.
- Add 2.5 ml of 10X PBS and 10 ml of additional DH₂O water.
- Check pH – should be around 7.0 – 7.5. If not start over.

**Make fresh every time (use w/in 24h, store at 4°C and bring to room temp before using)
Keep away from LIVE material, only work with under the hood!!**

2/28/2006

Elisha Wood-Charlson

Rodriguez-Lanetty, M., Krupp, D.A. and Weis, V.M. (2004) *Mar Ecol Prog Ser* 275: 97-102.

Procedure

A. Titer algae (see algal cell count protocol for more detail)

1. Use clean, recently extracted algae
2. If really dark brown, dilute to a pale brown color (keep track of dilutions)
3. Place coverslip on hemocytometer
4. Mix sample well
5. Load both sides
 - a. 10 μ l with a 10 μ l pipette into the hemocytometer
 - b. Load at a 45° angle w/ constant speed – should fill uniformly
6. Count squares depending on density
 - a. Dense – count middle only
 - b. Medium – count corners and middle
 - c. Light – count all nine
7. Record # algae counted in # of squares
8. Repeat with other side
9. Repeat load and counts one more time
10. Calculate cells/ml
 - a. (# cells/# squares) x (# squares) x (10,000) = # cells/ml
 - b. Average the 4 replicates to get a final conc. of cells/ml
11. If diluted original solution, calculate cells/ml in original solution
12. Titer to 10⁶ cells/ml for infections (inoculate w/ 1ml – 10⁶ total algae)

ml	l	ml
l	dml	l
ml	l	ml

B. Infecting *Fungia scutaria* larvae

1. Day 3: examine larvae under microscope to determine if mouths are beginning to form.

Do not infect before Day 4

2. Morning: conc. larvae during cleaning (see protocol “Cleaning *Fungia scutaria* larvae”)
 - a. Choose a populated, healthy larval bowl
 - b. Clean, but do not add fsw to rinsed larvae
 - c. Choose another bowl, clean and add larvae to previous bowl
 - d. Repeat if necessary – should be thick!
 - e. Use within 4-6hrs
3. Afternoon: infect w/ titered algae
4. Thaw *Artemia* sp., homogenize and filter to create an ‘essence of *Artemia*’
5. Add to each well in 6-well culture dish
 - a. 4ml fsw
 - b. 5ml conc. larvae
 - c. 1ml algae (= 10⁶ cells minimum) + ‘essence of *Artemia*’
6. Cover and place in watertable on bricks so immersed but not submerged
7. Feed for 3h

C. Rinsing infected *Fungia scutaria* larvae

1. After 3h, remove 6-well plates from watertable
2. Set up a clean 6-well plate, fill 4 wells w/ fsw
3. Using a clean plastic pipette, withdraw larvae from well and squirt onto 50 μ m cell filter over empty well

2/28/2006

Elisha Wood-Charlson

Rodriguez-Lanetty, M., Krupp, D.A. and Weis, V.M. (2004) *Mar Ecol Prog Ser* 275: 97-102.

4. Rinse larvae on the filter by dipping and squirting w/ fsw
5. Rinse original 6-well plate with di-water
6. Replace larvae to original well by inverting filter and squirting w/ fsw
7. Cover, return to watertable and repeat
8. When finished – this is TIME 0

D. Collecting and fixing *Fungia scutaria* larvae

1. At set time point, remove 1ml of larvae with plastic pipette and place in a LABELED 1.5ml tube
2. Spin 1min, 1000xg
3. Remove sup to conc. larvae
4. Fill with paraformaldehyde
5. Parafilm and store at 4°C

Comments: